

# PDG Fuze<sup>®</sup> M50

Manual tilt-in-space wheelchair



permobil

# Full tilt with optimal wheel access

Experience up to 50° of tilt for effective pressure relief and positioning, all while maintaining ideal drive wheel access for confident, independent manoeuvring.

Tackle everyday obstacles with ease thanks to new Tip Activation, which lifts the rear castors to provide up to 4" of front castor clearance, enhancing access in challenging environments.

**Tip activation:** Up to 4" of front castor clearance

**Tilt mechanism:** Up to 50° of tilt

**Self-propel:** Wheel position allows optimal access for hand propulsion

**6-wheel:** Mid-wheel drive approach for tight turning radius



## Specifications

|                                |              |
|--------------------------------|--------------|
| <b>Dynamic tilt range</b>      | Up to 50°    |
| <b>Seat width</b>              | 15"-24"      |
| <b>Seat depth</b>              | 15"-22"      |
| <b>Seat to floor height</b>    | 13"-19"      |
| <b>Overall width</b>           | + 8.5"       |
| <b>Transport weight</b>        | 24 kg        |
| <b>Maximum user weight</b>     | 181 kg       |
| <b>Frame material</b>          | Aluminium    |
| <b>Power tilt option</b>       | No           |
| <b>Take-Apart</b>              | No           |
| <b>Hand or foot propulsion</b> | Hand         |
| <b>Transit tiedowns</b>        | Up to 136 kg |

## Frame colours


### High gloss paint




Red pepper




Neon pink




High gloss black




Seafoam green



Lavender dusk




Blue mist



Daffodil yellow

### Pearlescent paint




Midnight blue pearl



White pearl


### Matte textured paint



Flame orange




OD green



Hyper blue

### Matte paint




Titanium grey




Black lava

### Candy paint




Pomegranate red




Indigo purple


### Metallic paint




Tangerine metallic




Acid green Metallic



Ocean blue metallic



Black Galaxy



Silver Streak