

PDG Mobility Bentley LT-R

The Bentley LT-R offers a reversed wheel layout for improved access to drive wheel, easier self-propelling and superior manoeuvrability, even when tilted.

The Bentley LT-R's front pivot design maintains the front seat-to-floor height and can allow users to keep their feet on the group even when the chair is fully tilted – a perfect solution for foot propelling or standing transfers.



Exclusive features

- Reverse wheel layout for improved independent self-propulsion
- Ability to hand and foot propel in a tilted position with no knee-rise
- User tilt activation for additional independence
- Lightweight wheelbase for good manoeuvrability

Available options (see order form)

- Transit tie downs WC19 & ISO7176-19
- 30° manual dynamic recline
- Power tilt
- Easy-Out front rigging
- Drum brakes



Hand Propulsion



Foot Propulsion



Transit Tie Downs



Self-Tilt



Desk/Table Accessible

Specifications

Dynamic tilt range	Up to 20°
Seat width	14" - 22"
Seat depth	16" - 22"
Seat to floor height	12" - 19"
Overall width	+10.5"
Transport weight	15 kg
Maximum user weight	113 kg
Frame material	Aluminium
Power tilt option	Yes
Take-Apart frame	No
Transit tie downs	Up to 113 kg

permobil

Scan the QR code to learn more about PDG Mobility.

