

# PDG Mobility Bentley

The Bentley offers up to 20° tilt for comfort and postural support. Its steel frame provides a high level of stability and durability.

The Bentley's front pivot design maintains the front seat-to-floor height and can allow users to keep their feet on the ground even when the chair is fully tilted – a perfect solution for foot propelling or standing transfers.



## Exclusive features

- Ability to hand and foot propelling in a tilted position with no knee-rise
- User tilt activation for additional independence
- Short wheelbase for manoeuvrability
- Adjustable rear axle position for centre of gravity changes

## Available options (see order form)

- Transit tie downs WC-19 & ISO-7176-19
- 30° manual dynamic recline
- Power tilt
- Drum brakes
- Floor jacks
- Easy-Out front rigging



Hand Propulsion



Foot Propulsion



Transit Tie Downs



Self-Tilt



Desk/Table  
Accessible

## Specifications

Dynamic tilt range	Up to 20°
Seat width	14" - 32"
Seat depth	16" - 22"
Seat to floor height	13" - 20"
Overall width	+9.5" - 10.5"
Transport weight	23 kg
Maximum user weight	204 kg
Frame material	Steel
Power tilt option	Yes, up to 113 kg
Take-Apart frame	No
Transit tie downs	Up to 136 kg

permobil

Scan the QR code to learn  
more about PDG Mobility.

