

# PDG Mobility Fuze JR

The Fuze JR is a paediatric 50° tilt chair featuring a unique centre of gravity pivot tilt system and adjustable seat frame as standard.



## Exclusive features

- A unique centre of gravity pivot tilt design
- Effortless and maintenance-free tilt operation
- Compact frame designed specifically for paediatrics
- Adjustable seat frame as standard

## Available options (see order form)

- Transit tie downs WC-19 & ISO-7176-19
- 30° manual dynamic recline
- Take-Apart frame
- Drum brakes

## Specifications

Dynamic tilt range	0° to 50°
Seat width	12" - 18"
Seat depth	13" - 18"
Seat to floor height	13" - 20"
Overall width	+7.5" - 8.5"
Transport weight	19 kg
Maximum user weight	113 kg
Frame material	Aluminium
Power tilt option	No
Take-Apart frame	Yes
Transit tie downs	Up to 113 kg



Hand Propulsion



Transit Tie Downs

**permobil**

Scan the QR code to learn more about PDG Mobility.

