

PDG Mobility Fuze M50

The M50 provides the benefits of 50° positioning and pressure-relief tilt plus optimal access for independent self-propelling.

Obstacles can be climbed easily due to the new Tip Activation feature that allows for up to 4" of front castor clearance.



Exclusive features

- A unique 6-wheel design with the main drive wheel in the middle for optimal self-propelling access and a very tight turning circle
- Effortless and maintenance-free pivot tilt operation

Available options (see order form)

- Transit tie downs WC-19 & ISO-7176-19
- 30° manual dynamic recline
- Adjustable width frame



Hand Propulsion



Transit Tie Downs

Specifications

| | |
|----------------------|--------------|
| Dynamic tilt range | 0° to 50° |
| Seat width | 15" - 24" |
| Seat depth | 15" - 22" |
| Seat to floor height | 13" - 19" |
| Overall width | +8.5" |
| Transport weight | 24 kg |
| Maximum user weight | 181 kg |
| Frame material | Aluminium |
| Power tilt option | No |
| Take-Apart frame | No |
| Transit tie downs | Up to 136 kg |

permobil

Scan the QR code to learn more about PDG Mobility.

