


PDG Mobility Stellar HD

The Stellar HD is a bariatric tilt wheelchair offering up to 45° of tilt and 272 kg weight capacity. It features a heavy-gauge steel reinforced frame and multiple gas struts to ensure effortless tilt operation.



Available options (see order form)

- 30° manual dynamic recline
- Dynamic rocker back
- Easy-out front rigging



Specifications

Dynamic tilt range	0° to 45° (25° at 272 kg)
Seat width	22" - 32"
Seat depth	18" - 22"
Seat to floor height	16" - 18"
Overall width	+11"
Transport weight	28 kg
Maximum user weight	272 kg
Frame material	Steel
Power tilt option	No
Take-Apart frame	No
Transit tie downs	No



Hand Propulsion



Desk/Table
Accessible

permobil

Scan the QR code to learn
more about PDG Mobility.

