

# PDG Mobility Stellar Impact

The Stellar Impact is designed for the most demanding wheelchair users.

Offering up to 45° tilt, its ultra-tough steel construction provides maximum durability for high tone, highly agitated, and/or continuously rocking occupants.



## Available options (see order form)

- Transit tie downs WC19 & ISO7176-19
- 30° manual dynamic recline
- Dynamic rocker back
- Drum brakes
- Floor jacks



## Specifications

Dynamic tilt range	0° to 45°
Seat width	14" - 22"
Seat depth	17" - 22"
Seat to floor height	14" - 21"
Overall width	+9.5 - 10.5"
Transport weight	29 kg
Maximum user weight	113 kg
Frame material	Steel
Power tilt option	No
Take-Apart frame	No
Transit tie downs	Yes



Hand Propulsion



Transit Tie Downs



Desk/Table  
Accessible

permobil

Scan the QR code to learn  
more about PDG Mobility.

