

Manual wheelchairs

TiLite CR1

permobil



Key features

- Exceptionally lightweight frame
- ± 1 " Centre of Gravity fine tuning adjustment
- Integrated, fold-down backrest for easy transport
- Individually made, just for you

Custom profiled sideguard

- Supports backrest without need for backrest hinges
- Individually designed for tailored fit

[12.0096, 12.0116]



Performance components

- Specially designed castor bearings for optimal manoeuvrability and weight savings
- PBO fibre and titanium reinforcements
- New carbon fibre castor wheels

Dual-angle front end


- Designed to increase environmental access
- Optimised castor placement





Component Colours

Anodised Options



Red

Orange

Blue

Purple

Green

Gold

Charcoal
Grey

COG Bracket & Trim Rings



Charcoal
Grey

Brushed
Aluminium

Specifications

Transport weight	3.9 kg	Camber	0°, 2°, 4°
Seat width	13" – 18"	Centre of gravity adjustment	2"
Seat depth	13" – 19"	Backrest angle	80° - 101°
Front seat height	17.5" – 21"	Compatible with SmartDrive?	Yes
Seat slope	7°, 10°, 13°	Weight limit	105 kg
Front frame angle	85°	Warranty	Limited
Seat back height	10" – 17"		