


Power assist

Progeo Mototronik



permobil




Mototronik is a stylish and intuitive front-mounted power assist device designed to improve your access and independence.

- Cover longer distances with a range of up to 50km*
- Cruise at higher speeds of up to 20km/h
- Front castors are elevated to make it easier to ride over challenging surfaces
- Easy on-screen profile and speed adjustment
- On-board programming to adjust your preferred settings if necessary
- Cruise control or throttle at your desired speed

*Dependent on configuration

Fully programmable settings with easy speed adjustment and cruise control



LEARN MORE

Optional LED headlight and detachable battery for easy charging

Available with Standard or Quadriplegic handlebars

12" or 14" wheel options



Removable stand and quick, intuitive attachment system

Compatible with most rigid and folding manual wheelchairs

Specifications

Power (12" wheel)	250w (max 750w)
Power (14" wheel)	1000w (max 1500w)
Motor voltage (12" & 14")	48V
Available speeds	up to 20km/h
Power Plus battery (11.6Ah)	up to 50km*
Standard battery (5.8Ah)	up to 25km*
Airline approved battery (FLY2.9Ah)	up to 15km*

Maximum towing load	120kg**
Maximum rated slope	10%
Transport weight (12" wheel)	from 10kg
Transport weight (14" wheel)	from 14kg
Expected service life	5 years or 25,000km

Optional extras include:
 Grip throttle, handlebar extension, LED headlight, rear-view side mirrors, travel battery, removable basket, removable weights, and travel bag.

* With a charged battery, moderate speed, flat course and user weight 75kg
 ** User's weight plus wheelchair weight