


Transportation Made Easy

The Take-Apart frame allows attendants to take the Fuze™ T50, Fuze T50n or Fuze JR apart in seconds without tools. This reduces the weight of the heaviest piece of the wheelchair to approximately 13kg, perfect for convenient transportation in the boot of a car or SUV.



Our innovative chair designs deliver an enhanced quality of life for both users and caregivers. Long product life, adjustability, and durability create value for all stakeholders, and our expertise in positioning wheelchairs ensures an excellent fit every time, no matter how basic or complex.

Seating and Positioning

Permobil provides a wide range of seating and positioning products to address the unique needs of individuals using PDG wheelchairs.

Comfort and ROHO® offer a variety of cushions, back supports and accessories that are designed for pressure redistribution, stability, and comfort that can be tailored to meet the individual needs of users.



PDG Mobility™

Comparison matrix



permobil

PDG Manual Wheelchairs



Fuze™ T50

Fuze T50n

Fuze T50 Junior

Fuze M50

Stellar GL

Stellar GLT

Stellar Impact

Stellar HD

Bentley™

Bentley LT

Bentley LT-R

Eclipse

	Fuze™ T50	Fuze T50n	Fuze T50 Junior	Fuze M50	Stellar GL	Stellar GLT	Stellar Impact	Stellar HD	Bentley™	Bentley LT	Bentley LT-R	Eclipse
Dynamic tilt range	0° to 50°	0° to 50°	0° to 50°	0° to 50°	-5° to 45°	-5° to 45°	0° to 45°	0° to 45° (25° at 272kg)	0° to 20°	0° to 20°	0° to 20°	None
Seat width	12" - 24"	16" to 24"	12" - 18"	15" - 24"	14" - 32"	14" - 22"	14" - 22"	22" - 32"	14" - 32"	14" - 22"	14" - 22"	20" - 48"
Seat depth	15" - 22"	15" - 22"	13" - 18"	15" - 22"	16" - 22"	16" - 22"	17" - 22"	18" - 22"	16" - 22"	16" - 22"	16" - 22"	16" - 24"
Seat to floor height	13" - 20"	16" - 20"	13" - 20"	13" - 19"	13" - 20"	13" - 20"	14" - 21"	16" - 18"	13" - 20"	13" - 20"	12" - 19"	14" - 20"
Overall width*	+7.5"	+4.0"	+7.5"	+8.5"	+9.5"	+9.5"	+9.5"	+11"	+9.5"	+9.5"	+10.5"	+7.5"
Transport weight	20kg	20kg	19kg	24kg	23kg	15kg	29kg	28kg	23kg	15kg	15kg	20kg
Maximum user weight	181kg	181kg	113kg	181kg	204kg	113kg	113kg	272kg	204kg	113kg	113kg	454kg
Frame material	Aluminium	Aluminium	Aluminium	Aluminium	Steel	Aluminium	Steel	Steel	Steel	Aluminium	Aluminium	Steel
Power tilt option	Yes, up to 159kg	Yes, up to 159kg	Yes	No	Yes, up to 113kg	Yes	No	No	Yes, up to 113kg	Yes	Yes	No
Take-Apart frame option	Yes	Yes	Yes	No	No	No	No	No	No	No	No	No
Transit tiedowns	Up to 136kg	Up to 136kg	Up to 113kg	Up to 136kg	Up to 136kg	Yes	Yes	No	Up to 136kg	Up to 113kg	Up to 113kg	Up to 272kg

*Overall width measurement is in addition to the seat width.