

Permobil

# Progeo Duke

Individualized performance with a choice of carbon



Carbon



Rigid



Minimized  
weight



Adjustable

permobil

# Individualized with the room for adjustments

Progeo Duke is a rigid carbon wheelchair with a fixed front frame that offers seating and drive performance adjustability. The fixed front and the carbon material enhance the maneuverability resulting in a more durable and agile ride. This makes the Progeo Duke a solution for users favoring the benefits of carbon.



## Progeo Duke

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
Seat height	Front from 44,5 to 53 cm (adjustable) Rear from 38 to 46 cm (adjustable)
Footplate distance	From 33 to 49 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	90°
Total length	85 cm approx. (standard configuration)
Total width	Seat width + 17 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 5,2 kg approx.
Wheelchair weight	Starting from 7,8 kg approx.
Maximum load	125 kg

Permobil reserves the right to change specifications and design without prior notice. \*More color options available.

Personalise your ride\*

MATTE  
BLACK

BRIGHT  
LIME

BRIGHT  
FUCHSIA

BRIGHT  
BLUE

BRIGHT  
RED

Actual product colors may vary.