

Permobil

Progeo Ego

Designed for independence



Aluminium



Folding



Adjustable

permobil

The performance of a rigid chair in a folding frame

Progeo Ego is designed around a horizontal folding system that locks. The locking system makes it perform like a rigid chair. Combined with adjustability of the seating and the drive performance the user will have a well balance and energy efficient ride. Additionally the open frame design and folding backrest make it truly compact for ease of transportation.



Progeo Ego

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
Seat height	Front from 40 to 58 cm (adjustable) Rear from 38 to 45 cm (adjustable)
Footplate distance	From 35 to 54 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	Aluminium 85° - 90° - 100° Carbon 85° - 90° - 95°
Total length	88 cm approx. (standard configuration)
Total width	Seat width + 17 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2°
Weight without rear wheels	Starting from 7,2 kg approx.
Wheelchair weight	Starting from 9,8 kg approx.
Maximum load	125 kg

Permobil reserves the right to change specifications and design without prior notice. *More color options available

Personalise your ride*

POLISHED
BLACK

BRIGHT
LIME

BRIGHT
RED

BRIGHT
BLUE

BRIGHT
SILVER

Actual product colors may vary.