

Permobil

Progeo Exelle Vario

Everyday chair designed with the user in mind



Aluminium



Foldable



Adjustable

permobil

Compact and elegant with numerous adjustments

Progeo Exelle Vario is a lightweight, folding wheelchair designed for the active user. It's contoured folding frame and the swing-in/out footrests enable easy transfers and storage. A variety of options to support seating and positioning such as the three legrest frame angles can address knee flexion/extension and spasticity management. The adaptability built into the frame makes this a versatile solution for many different users.



Progeo Exelle Vario

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
Seat height	Front from 41 to 56 cm (adjustable) Rear from 36 to 53 cm (adjustable)
Footplate distance	From 32 to 51 cm (adjustable)
Backrest height	From 30 to 54 cm (adjustable)
Front frame angle	95° - 110° - 120°
Total length	90 cm approx. (std. config.) 70 cm approx. (without footplates) (std. config.)
Total width	Seat width + 17 cm (with 0° camber)
Encumbrance	26 cm approx. (with 0° camber)
Backrest angle	90° to the seat (with standard rear frame) Adjustable from 81° to 99° to the seat (with Dynamic rear frame)
Camber	0° - 1,5° - 3° (with standard rear frame) 0° - 1,5° - 2° - 4° (with Dynamic rear frame)
Weight without rear wheels	Starting from 8,0 kg approx.
Wheelchair weight	Starting from 10,6 kg approx.
Maximum load	125 kg

Permobil reserves the right to change specifications and design without prior notice. *More colour options available.

Personalise your ride*

- BLACK
- BRIGHT LIME
- BRIGHT FUCHSIA
- BRIGHT BLUE
- BRIGHT SILVER

Actual product colours may vary.