

Permobil

Progeo Joker

Adjustability meets performance



Aluminum



Rigid



Adjustable

permobil

On demand adjustability

Progeo Joker is combining the benefits of a rigid frame with adjustability of the seating and the drive performance. This gives the user a balanced and responsive ride. The adjustability makes this a preferred choice for any user with changing needs.



Progeo Joker

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
Seat height	Front from 40 to 59 cm (adjustable) Rear from 38 to 46 cm (adjustable)
Footplate distance	From 35 to 54 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	Aluminium 85° - 90° - 100° Carbon 85° - 90° - 95°
Total length	88 cm approx. (standard configuration)
Total width	Seat width + 17 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 5,4 kg approx.
Wheelchair weight	Starting from 8,0 kg approx.
Maximum load	125 kg

Personalise your ride*

BLACK

BRIGHT LIME

BRIGHT FUCHSIA

BRIGHT BLUE

BRIGHT RED

Actual product colors may vary.

Permobil reserves the right to change specifications and design without prior notice. *More color options available