

Permobil

# Progeo Joker Carbon

Ultra-light Carbon Fibre Rigid Wheelchair



Carbon



Rigid



Adjustable

permobil

## A carbon frame that handles more

JOKER CARBON, a new evolution of the model Joker characterised by its monocoque carbon fibre frame that remarkably reduces the weight while improving manoeuvrability. Weight in using configuration starting from 7.8 kg.



## Progeo Joker Carbon

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
Seat height	Front from 45 to 53 cm (adjustable) Rear from 36 to 48 cm (adjustable)
Footplate distance	From 35 to 49 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	Carbon 90°
Total length	88 cm approx. (standard configuration)
Total width	Seat width + 17 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 5,0 kg approx.
Wheelchair weight	Starting from 7,8 kg approx.
Maximum load	125 kg

Personalise your ride

BLACK

BRIGHT ORANGE

BRIGHT RED

BRIGHT FUCHSIA

BRIGHT BLUE

Actual product colours may vary

Permobil reserves the right to change specifications and design without prior notice. \*More colour options available.