

Permobil

Progeo Noir 2.0

Performance and design in perfect harmony



Carbon



Rigid



Minimized
weight



Made-to-
measure

permobil

The perfect fit

Progeo Noir 2.0 combines the benefits of carbon with the superior manoeuvrability of a made to measure wheelchair. This together gives the user an optimal balanced ride allowing for energy saving throughout the day.



Progeo Noir 2.0

| | |
|----------------------------|---|
| Seat width | 33 - 36 - 39 - 42 - 45 cm |
| Seat depth | 35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm |
| Seat height | Made-to-measure |
| Footplate distance | Made-to-measure (standard) Adjustable (on request) |
| Backrest height | From 24,5 to 47 cm (adjustable) |
| Front frame angle | 85° - 90° - 95° |
| Total length | 85 cm approx. (standard configuration) |
| Total width | Seat width + 17 cm (with 0° camber) |
| Backrest angle | From 74° to 94° to the seat (adjustable) |
| Camber | 0° - 2° - 4° |
| Weight without rear wheels | Starting from 4,0 kg approx. |
| Wheelchair weight | Starting from 6,6 kg approx. |
| Maximum load | 125 kg |

Personalise your ride

CARBON
ESSENCE

LIME
ESSENCE

RED
ESSENCE

BLUE
ESSENCE

OPAQUE
ESSENCE

Actual product colors may vary.

Permobil reserves the right to change specifications and design without prior notice.