

Permobil

Progeo Tekna Tilt 2.0

Individualized positioning meets mobility



Aluminium



Foldable



Adjustable



Tilt

permobil

Variable seating and positioning with compact folding

Progeo Tekna Tilt 2.0 is a compact, folding wheelchair offering tilt, optional recline and elevating legrest to accommodate a broad range of seating and positioning needs. It is one of the very few tilt-in-space wheelchairs that folds into a compact set-up for transportation or storage. Users can reposition according to their needs and select either an active and passive mobility version.



Progeo Tekna Tilt 2.0

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
Seat height	Min. 42 Max 48 cm. According to front and rear wheels. See order form
Footplate distance	From 27 to 45 (adjustable)
Backrest height	From 35 to 55 cm (adjustable)
Front frame angle	Standard frame 95° - 110° - 120° Frame with elevating leg rest from 100° to 180°
Total length	90 cm approx. 70 cm approx. without footplates (std. config.)
Total width	Seat width + 22 cm (with 300 mm rear wheels)
Encumbrance	32 cm (with 300 mm rear wheels)
Backrest angle	From 78° to 94° to the seat (adjustable)
Tilt-in-space	From 0° to 45°
Weight without rear wheels	Starting from 12 kg approx. (without footplates, rear wheels, armrests and anti-tippers)
Wheelchair weight	Starting from 17,5 kg approx.
Maximum load	125 kg

Permobil reserves the right to change specifications and design without prior notice. *More color options available.

Personalise your ride*

BLACK

BRIGHT LIME

BRIGHT FUCHSIA

BRIGHT BLUE

BRIGHT SILVER

Actual product colors may vary.