

permobil



Ready Set  
Roll!

***K300***  
PS Jr.



### Powerful Suspension

Provides a smoother ride, even over rough terrain, to help you handle the adventures ahead.



### Highly Configurable

Allows you to maintain proper positioning as you grow.



### Multi-Motion Armrest

Perfectly positioned armrest with unique rotation and flip back movements that allow for easier transfers.

### Colour Options



Volcano Red



Pop Star Pink



Cobalt Blue



Purple Rain



Phantom Gray



Onyx Black

### Additional Hub Colour Options



Chili Pepper Red



Titan Blue



Extreme Yellow



Vintage Wine



Crystal Blue

### Specifications

|                            |  |
|----------------------------|--|
| Maximum User Weight        | 75 kg  |
| Maximum Speed              | 8 km/h or 10 km/h                                  |
| Range                      | 25 km  |
| Base Width                 | 24"  |
| Base Length                | 40"  |
| Minimum Turning Radius     | 25"  |
| Ground Clearance           | 3"   |
| Recommended Battery Type   | Group 34   |
| Weight Including Batteries | 123.3 kg   |
| Drive Electronics          | 90A R-Net  |
| Maximum Slope              | 10° up, 6° down <sup>1</sup>                       |
| Seat-to-Floor Height       | 16.5" (elevator)<br>16.5", 17.25", 18" (seat tube) |

|                        |   |
|------------------------|---|
| Optional Seat Elevator | 8"  |
| Tilt Options           | 0° – 45° (power)<br>0° – 45° (manual)                   |
| Manual Recline         | 85° – 120°  |
| Legrest Elevation      | 70° (swing-away)<br>90° – 150° (center mount elevating) |
| Backrest Height        | 16.5" – 22.5"   |
| Seat Depth             | 10" – 18"   |
| Seat Width             | 11" – 16"   |
| Armrest Pad Lengths    | 10", 13"  |

<sup>\*</sup>Actual driving range will vary based on driving and battery conditions.  
<sup>1</sup>Refer to Owner Manual for safe driving practices.