

permobil

Elevate play

M300
PS JR.





Highly adjustable seat that grows as fast as you do.



Catlink Suspension

Provides a smoother ride, even over rough terrain, to help you handle the adventures ahead.



Highly Configurable


Allows you to maintain proper positioning as you grow.




Multi-motion Armrest

Perfectly positioned armrest with unique rotation and flip back movements that allow for easier transfers.

Shroud Colour Options



Additional Hub Colour Options



Specifications

Maximum User Weight	75 kg
Maximum Speed	10 km/h
Range	20 km*
Base Width	24"
Base Length	36.5"
Minimum Turning Radius	20"
Ground Clearance	3"
Recommended Battery Type	Group 24, Group 34
Weight Including Batteries	142 kg G34 / 151 kg G24
Drive Electronics	120A R-Net
Maximum Slope	9° up, 9° down ¹
Seat-to-Floor Height	17" (elevator) 17", 17.75", 18.5" (seat tube)

Optional Seat Elevator	8"
Tilt Options	0° – 45° (power) 0° – 45° (manual)
Manual Recline	85° – 120°
Legrest Elevation	70° (swing-away) 90° – 150° (centre mount elevating)
Backrest Height	16.5" – 22.5"
Seat Depth	10" – 18"
Seat Width	11" – 16"
Armrest Pad Lengths	10", 13"

*Actual driving range will vary based on driving and battery conditions.
¹Refer to Owner Manual for safe driving practices.